

Make informed decisions quickly and conveniently.

Alere Oratect®

Oral Fluid Drug Screen Device



The first U.S. in vitro diagnostic oral fluid drug test device.


Knowing now matters.™

Alere Oratect® Oral Fluid Drug Screen Device


Features and Benefits

- Rapidly screen 6 illicit drugs and prescription drugs in oral fluid
- Simple procedure: collect saliva and test in 1 step—results in five minutes
- Convenient and noninvasive rapid screening
- Eliminates need for gender specific observation
- Reduces sample adulteration or contamination
- Includes transport tubes for confirmation testing
- For in vitro diagnostic use in the U.S.


Test in 3 Easy Steps



1
Remove plastic cap at bottom. Collect saliva sample.



2
Gently swab the inside of the mouth and tongue with pad.



3
Replace cap on device. Lay flat. Read test results at 5 minutes.



PRODUCT PROCEDURE NOTE: Refer to product insert for complete instructions, limitations, precautions, and warnings. For in vitro diagnostic use and it can be used at the point of care site. This device is a one-step, lateral-flow immunoassay device and provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography-mass spectrometry (GC-MS) or liquid chromatography-tandem mass spectrometry (LC-MS/MS) are the preferred confirmatory methods. The samples for confirmatory testing should be collected with the oral fluid confirmation tube provided.

Drug screen options and cutoff levels

Test	Calibrator	Cutoff
Amphetamine (AM)	d-Amphetamine	50 ng/mL
Cocaine (CO)	Cocaine	20 ng/mL
Marijuana (TH)	Δ -9-Tetrahydrocannabinol	40 ng/mL
Methamphetamine (ME)	d-Methamphetamine	50 ng/mL
Opiate (OP)	Morphine	40 ng/mL
Phencyclidine (PC)	Phencyclidine	10 ng/mL

ITEM #: AOT-6

25 devices per box. Includes 2 collection tubes for confirmation testing.

