

Dependable results and simple to administer.

iCassette[®]

Urine Drug Screening Device.




Substance abuse testing with more **substance**.

How it works.

PRODUCT TRAINING AND
CERTIFICATION AVAILABLE ONLINE.

Features & Benefits

- » Rapidly screen up to 11 illicit and prescription drugs
- » FDA 510(k) Cleared; CLIA Waived 6 drug test available for clinical use
- » Simple procedure: collect and pipette to test—results in as little as 5 minutes
- » Test device remains in control of administrator
- » Long shelf-life



PRODUCT PROCEDURE NOTE: Refer to product insert for complete instructions, limitations, and warnings.

For in vitro diagnostic use only. This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgement should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

Urine Screening Device Cutoff Levels

Test	Calibrator	Cutoff
Amphetamine (AMP)	d-Amphetamine	1,000 ng/mL
Barbiturates (BAR)	Secobarbital	300 ng/mL
Benzodiazepines (BZO)	Oxazepam	300 ng/mL
Buprenorphine (BUP)	Buprenorphine	10 ng/mL
Cocaine (COC)	Benzoylcegonine	300 ng/mL
Marijuana (THC)	11-nor- Δ^9 -THC-9-COOH	50 ng/mL
Methadone (MTD)	Methadone	300 ng/mL
Methamphetamine (mAMP)	d-Methamphetamine	1,000 ng/mL
Methylenedioxymethamphetamine (MDMA) Ecstasy	d,l-Methylenedioxymethamphetamine	500 ng/mL
Opiate (MOP)	Morphine	300 ng/mL
Opiate (OPI)	Morphine	2,000 ng/mL
Oxycodone (OXY)	Oxycodone	100 ng/mL
Phencyclidine (PCP)	Phencyclidine	25 ng/mL
Propoxyphene (PPX)	Propoxyphene	300 ng/mL
Tricyclic Antidepressants (TCA)	Nortriptyline	1,000 ng/mL

