

Multiple configurations to meet your unique drug test needs.

iScreen[®] Dx Dip Card

+

Urine Drug Screening Device.



Featuring

Buprenorphine test and lower
Opiate cutoff levels (300 ng/mL)

*Made in the U.S.A.
For in vitro diagnostic use only.*


Make smart and informed decisions.

How it works.

PRODUCT TRAINING AND
CERTIFICATION AVAILABLE ONLINE.

Features & Benefits

- » CLIA-waived
- » Rapidly screen for a broad range of illicit and prescription drugs
- » 1-12 configurations available for a fully customizable urine testing program
- » Simple procedure: collect, dip and read—results in as little as five minutes
- » Test device remains in control of administrator



IMPORTANT PROCEDURE NOTE: Refer to product insert for complete instructions, limitations, and warnings. For in vitro diagnostic use only. Urine drug screen assay's provide only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) or liquid chromatography/tandem mass spectrometry (LC/MS/MS) are the preferred confirmatory methods. Clinical consideration and professional judgement should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

Single Drug Configurations (25 per box)

| Item # | iScreen® Dx (Dip Card) |
|------------|---|
| IS1 AMP-DX | Amphetamine (AMP) Single Dip Device |
| IS1 BAR-DX | Barbiturate (BAR) Single Dip Device |
| IS1 BUP-DX | Buprenorphine (BUP) Single Dip Device |
| IS1 BZO-DX | Benzodiazepine (BZO) Single Dip Device |
| IS1 COC-DX | Cocaine (COC) Single Dip Device |
| IS1 MET-DX | Methamphetamine (MET) Single Dip Device |
| IS1 MTD-DX | Methadone (MTD) Single Dip Device |
| IS1 MOP-DX | Morphine (MOP-300) Single Dip Device |
| IS1 OXY-DX | Oxycodone (OXY) Single Dip Device |
| IS1 PCP-DX | Phencyclidine (PCP) Single Dip Device |
| IS1 MDM-DX | Ecstasy (MDMA) Single Dip Device |
| IS1 TCA-DX | Tricyclic Antidepressants (TCA) Single Dip Device |
| IS1 THC-DX | Marijuana (THC) Single Dip Device |

Multi-drug Configurations (25 per box)

| Item # | iScreen® Dx (Dip Card) |
|---------------|---|
| IS5 M-DX | AMP, COC, MET, MOP-300, THC |
| IS10-DX | AMP, BAR, BZO, COC, MET, MOP-300, MTD, OXY, TCA, THC |
| IS10-22015-DX | AMP, BAR, BZO, COC, MET, MDMA, MOP-300, MTD, PCP, THC |
| IS10BUP-DX | AMP, BAR, BUP, BZO, COC, THC, MET, MOP-300, MTD, OXY |
| IS12-DX | AMP, BAR, BZO, COC, MET, MDMA, MOP-300, MTD, OXY, PCP, TCA, THC |

